


- **Halogen-Metal exchange**

The Halogen-metal exchange reaction was discovered in the late 1930s by Gilman⁵ and Wittig⁶. The reaction is often used to prepare vinyl- and aryllithium compounds from the more reactive alkyllithium species⁷.

- **Transmetallation**

The organolithium compounds are very often used to prepare other metallorganic compounds through the transmetallation reaction:



This reaction has been used to prepare the important organo-copper^{8, 9, 10} (A) and organo-titanium¹¹ (B) and many other metallorganic compounds¹² which have often higher selectivity than the organolithium compounds.

- **Anionic Polymerisation**

A major industrial use of alkyllithium compounds, specifically n-Butyllithium, is the catalysis of the anionic polymerization of butadiene¹³, isoprene and styrene.

In reactions involving organolithium reagents it is essential to exclude water and air, which rapidly destroy the reagent. Therefore, we provide these sensitive reagents in the industry-leading AcroSeal[®] packaging, preventing the degradation of the reagents.



Literature references:

1. Addition to Tosylhydrazione: A. G. Myers and M. Movassaghi, J. Am. Chem. Soc. 1998, 120, 8891-8892.
2. R.C.Larock, Comprehensive Organic Transformations VCH, New York 1989. p 512.
3. R.C.Larock, Comprehensive Organic Transformations VCH, New York 1989. p 792.
4. Leads with excess of RLi to ketones: T.M.Bare, H.O.House, Org.Synth. 49 (1969) 81.
5. H.Gilman, W.Langham, A.L.Jacoby, J.Am.Chem.Soc. 1939, 61, 106.
6. G.Wittig, U.Pockels, Chem.Ber. 1938, 71, 1903.
7. F.Leroux, M.Schlosser, Angew. Chem. 2002 114, 4447.
8. E. Erdik 1984, 40, 641; J.Lindley, Tetrahedron 1984, 40, 1433; B. H. Lipshutz, Sengupta, S. Org. React. 1992, 41, 135; Y. Yamamoto, Angew.Chem.Int.Ed. 1986, 25, 947; B. H. Lipshutz Synthesis 1987, 325; Marshall, J. A. Chem. Rev. 1989, 89, 1503; E.Nakamura, Synlett 1991, 539; Ibuka, T.; Yamamoto, Y. Synlett 1992, 769; Wipf, P. Synthesis 1993, 537; Krause, N.; Gerold, A. Angew. Chem. Int. Ed. Engl. 1997, 36, 187-204.
9. B.H.Lipshutz, Synthetic Procedures Involving Organocopper Reagents, Organometallics in Synthesis, M.Schlosser (Ed.) John Wiley & Sons 1994.
10. J.F.Normant, Synthesis 1972 63.
11. M.T.Reetz, Titanium in Organic Synthesis in Organometallics in Synthesis, M.Schlosser (Ed.) John Wiley & Sons 1994.
12. Manganese: Organic Syntheses, CV 9, 328; Organic Syntheses, Vol. 76, 239; Zinc: Organic Syntheses, CV 8, 430; Cerium: T.Imamoto, Pure&Appl.Chem. 1990 747; V. Nair, J.Mathew, Chem. Soc. Rev. 1997, 26, 127.
13. H.L.Hsieh, J.Polym.Sci. 1963 A3, 153.

Organolithium Compounds

AC377491000	Acros Organics	n-Butyllithium, 2.2M solution in cyclohexane, AcroSeal®	100ML	109-72-8
AC377498000			800ML	
AC213351000	Acros Organics	n-Butyllithium, 2.5M solution in hexanes, AcroSeal®	100ML	109-72-8
AC213358000			800ML	
AC181271000	Acros Organics	n-Butyllithium, 1.6M solution in hexanes, AcroSeal®	100ML	109-72-8
AC181275000			500ML	
AC181278000			800ML	
AC378931000	Acros Organics	n-Butyllithium, 2.7M solution in toluene, AcroSeal®	100ML	109-72-8
AC378938000			800ML	
AC187541000	Acros Organics	sec-Butyllithium, 1.3M sol. in cyclohexane/hexane (92/8), AcroSeal®	100ML	598-30-1
AC187548000			800ML	
AC396541000	Acros Organics	tert-Butyllithium, 1.9M solution in pentane, AcroSeal®	100ML	594-19-4
AC396548000			800ML	
AC445921000	Acros Organics	tert-Butyllithium, 2M (18%) solution in heptane, AcroSeal®	100ML	594-19-4
AC301651000	Acros Organics	n-Hexyllithium, 33 wt.% solution in n-hexane, AcroSeal®	100ML	21369-64-2
AC301658000			800ML	
AC377591000	Acros Organics	Isobutyllithium, 1.6M solution in heptane, AcroSeal®	100ML	920-36-5
AC181250500	Acros Organics	Lithium acetylide ethylenediamine complex, 85%	50GR	6867-30-7
AC181251000			100GR	
AC268831000	Acros Organics	Lithium diisopropylamide, 2M sol. in THF/n-heptane/ethylbenzene, AcroSeal®	100ML	4111-54-0
AC268838000			800ML	
AC432871000	Acros Organics	Lithium diisopropylamide mono(tetrahydrofuran), 1.5M solution in cyclohexane, AcroSeal®	100ML	123333-84-6
AC432878000			800ML	
AC380651000	Acros Organics	Lithium (trimethylsilyl)acetylide, 0.5M solution in THF/hexanes, AcroSeal®	100ML	54655-07-1
AC380658000			800ML	

AC188751000	Acros Organics	Methylithium, 1.6 M sol. in diethyl ether (\pm 5% w/v), AcroSeal®	100ML	917-54-4
AC188758000			800ML	
AC181291000	Acros Organics	Methylithium lithium bromide complex, 2.2M sol. in diethylether, AcroSeal®	100ML	332360-06-2
AC181298000			800ML	
AC445841000	Acros Organics	Methylithium, 3% solution in 2-MeTHF/cumene, AcroSeal®	100ML	917-54-4
AC445848000			800ML	
AC377451000	Acros Organics	(Trimethylsilyl)methylithium, 0.7M (10 wt%) solution in hexanes, AcroSeal®	100ML	1822-00-0

If you would like to know more about our extensive range of organolithium reagents, download the full product brochure from www.acros.com



To place an order, contact your local distributor.



©2015 Thermo Fisher Scientific Inc. All rights reserved.
All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.

Thermo Fisher Scientific
ENA 23, Zone 1, nr 1350
Janssen Pharmaceuticaaan 3a
2440 Geel
Belgium
www.acros.com